


Laser Safety Audits Of Research Labs



Johnny Jones

Laser-Professionals Inc.

Planning the Audit

- ✱ Write an objective statement.
- ✱ Create a data form that meets your needs.
- ✱ Ask all personnel about training.
- ✱ Make sketches of labs; take photos.
- ✱ Train audit team?

OPEN BEAM CONTROL MEASURES

ANSI Section 4.3.1.1

- ✱ Training
- ✱ Beam Control
- ✱ Laser Safety Eyewear
- ✱ Written Procedures
- ✱ Laser Controlled Area


First Question for the LSO: Have Training Requirements Been Met?

- ✱ Review Laser Safety Officer training.
- ✱ Review training records for laser personnel.
- ✱ Ask all personnel about training.

First Question In Every Lab: Where is the greatest hazard?

- ✱ Ignore entryway and interlocks in the beginning.
- ✱ Ask laser personnel to explain the setup, the hazards, and the controls.
- ✱ Review Written Procedures and SOPs.
- ✱ Evaluate Beam Control.
- ✱ Ask what could go wrong?

EXAMPLES OF GOOD BEAM CONTROL








**Surround hazardous
setups with barriers.**

**Block all reflections as near
their source as possible.**

CURBS ON OPTICAL TABLE




COMPUTERS IN RESEARCH LABS





Allowing a direct view from a computer workstation into a laser experimental setup increases the risk of eye exposure to reflected beams.




















TUNABLE ULTRASHORT PULSE SYSTEM




TI:SAPPHIRE ULTRASHORT PULSE SYSTEM



LASER SAFETY EYEWEAR




EYEWEAR LABELS




All eyewear must be labeled!

PLASTIC EYEWEAR CHARACTERISTICS



BEAM VIEWING DEVICES



Photos courtesy of



**WHO HAS PRIMARY RESPONSIBILITY
FOR LASER SAFETY ANY TIME A
CLASS 3B OR CLASS 4 LASER
IS OPERATED?**

**The person operating the laser
always has the primary
responsibility for all hazards
associated with laser use.**

Ask About Hazards During Alignment






Photo courtesy of Laser Safety Systems

SAFE BEAM ALIGNMENT

- Most beam injuries occur during alignment.
- Only trained personnel may align class 3B or class 4 lasers (**NO EXCEPTIONS!**)
- Laser safety eyewear is required for class 3B and class 4 beam alignment.
- ANSI **REQUIRES** approved, written alignment procedures for **ALL** class 4 laser alignment activities and recommends them for class 3B.

LASER PROTECTIVE BARRIERS





CLASS 4 ENTRYWAY CONTROLS

Section 4.3.10.2.2

1. *Non-Defeatable Entryway Controls*

- Doorway interlock is non-defeatable
- Training of authorized users only

2. *Defeatable Entryway Controls*

- Doorway interlock is defeatable
- Training of all personnel with access
- Barrier and eyewear at door

3. *Procedural Entryway Controls*

- No doorway interlock
- Training of all personnel with access
- Barrier and eyewear at door
- Visible or audible signal at doorway



No hazard at entryway

LABORATORY DOOR INTERLOCK



ENTRYWAY WARNING LIGHTS





DANGER



VISIBLE and/ or INVISIBLE LASER
RADIATION-AVOID EYE OR SKIN
EXPOSURE TO DIRECT OR
SCATTERED RADIATION.

ND:YAG 1064 nm

100 Watts Max. Average Power

Eyewear Required: OD ≥ 5 @ 1064 nm

CLASS 4 LASER

Controlled Area Warning Sign

Finish the Audit

- ✱ Send an audit report to all laser owners.
- ✱ Write a final audit report.
- ✱ Recommend improvements.

LASER PROFESSIONALS INC.



Experience Makes the Difference