



Virtual NMR at the HFMRF

David Hoyt, Ph.D.

Technical Group Leader/Senior Research Scientist

August 30, 2005











High Field Magnetic Resonance Facility

provides state-of-the-art nuclear magnetic resonance (NMR) and electron paramagnetic resonance (EPR) instrumentation for determining molecular structures that impact environmental remediation and biological health effects.

Instrumentation

- Twelve NMR spectrometers (ranging from 900 to 300 MHz) and one pulsed EPR spectrometer, with capabilities in high-field liquid-state, solid-state and micro-imaging techniques
- Combined optical and magnetic resonance microscope
- Low temperature probes for metallo-protein chemistry and structure
- High temperature probe for catalytic materials
- Virtual NMR capability to enable use and collaboration with EMSL scientists by remote users

Contact: David Hoyt, david.hoyt@pnl.gov (509) 373-9825









Science Areas

Solution-state

- **Structural Genomics**
- **Biomolecular complexes**
 - **Nucleic Acids/Protein**
 - **Protein-Protein complexes**
 - **Protein-Lipid complexes**

Solid-state

- Metalloproteins
- Solid-state materials research
 - •- catalytic, ceramics, amorphous materials.
 - •- emerging in biological solid-state research.
- Environmental chemistry
- Actinide chemistry

EPR

- Metalloproteins
- Membrane proteins and
- redox chemistry mechanisms
- Catalytic RNA's

MR Imaging/ Microscopy

- Imaging and metabolic "signatures" of cells/tissues and organisms.
- Flow and transport-
- Polymer aging
- Methods development
 - •- Slow-MAS
 - confocal microscopy





X EMSL

High Field Magnetic Resonance Facility

Proposal Process

The primary objective is to facilitate the best possible science.

- Open call twice a year (1800 recipients)
- Submitted online via the EMSL User System
- Proposals reviewed
 - EMSL staff
 - External panel
- Review criteria:
 - Scientific Merit
 - Appropriateness of the requested instrumentation
 - Relevance to the EMSL/DOE Mission
 - Contribution staff can have in bringing that science to fruition



Jan 15 April | July 15 October







Virtual NMR Access at HFMRF

- Secure, remote operation of spectrometers via internet
- Live consultations using collaboration tools
- Reduced travel costs and more flexible scheduling
- ~20% of facility users elect to use remote access tools



Publications:

Hoyt, D.W., et al. (2004) "Expanding your Laboratory by Accessing Collaboratory Resources". *Analytical and Bioanalytical Chemistry* **378 (1)**: 1408-1410.

Chin, G., et al. (2002) "Social Networks in the Virtual Science Laboratory". Communication of the Association for Computing Machinery 45 (8): 87-92.

Keating KA, et al. (2000) "Development and Use of a Virtual NMR Facility", *J. Mag. Res.*, <u>143</u> p. 172-183.







Virtual NMR Access at HFMRF

U. Washington - G. Varani Structure of telomerase RNA and telomeric proteins Montana State U. - Valerie Copie Three Dimensional Solution Structure of NoSL

Wayne State U. - J. Wang Structure of Apo-Lp III/HDL Particles

Cornell - Chih-hsin Cheng Properties of Surface Functional Groups of Black Carbon

Princeton - F. Hughson Intracellullar Protein Transport

CWRU - F. Sonnichsen Vanderbilt - Sanders 100 KDal Membrane-Bound Enzyme Complex

U. Central Florida - L. An Amorphous SiAICN Ceramics

EMSL

U.C. Davis -T. Dieckmann Structure Determination of Large RNA Oligomers

Burnham Institute -K. Ely Adaptor Proteins Implicated in Breast-Cancer Cell Drug

Resistance

Structure and Function of the Membrane Protein OFP16

Battelle

Arizona State U. - R. Nieman Los Alamos National Laboratories -Ryszard Michalczyk Studies of Ligand-Induced Conformational Changes in CD44 Receptor WWW EMSL.PNL.GOV

U. Kansas - C. Larive Wichita State - W. Carper Fulvic Acid-Heavy Metal Complexes

Pacific Northwest National Laboratory

Operated by Battelle for the U.S. Department of Energy



High Field Magnetic Resonance - Summary

- ► HFMRF: 12-NMR spectrometers, 1-pulsed EPR
- ▶ 900 MHz NMR

High Resolution liquids – structure of complexes
Solid State – optimal for quadrupolar nuclei and biological solids
63 mm bore – optimal for solid-state/ imaging projects

► Cold probes

2 - 600 MHz now available to users 800 MHz delivery (9/15/05)

- ► Next call for proposals (due 01/15/06)
- Remote access to operate instruments with staff



