


A Synchrotron X-ray National Facility for
Chemistry and Materials Research

NSF's ChemMatCARS: A Synchrotron X-ray National Facility for Chemistry and Materials Research




Introduction



Matthew Tirrell

*PI and Director of NSF's
ChemMatCARS, University of Chicago*

NSF's ChemMatCARS, Sector 15 at the Advanced Photon Source



Funded by NSF, located at the Advanced Photon Source (APS), Argonne National Laboratory, Illinois

Advanced Photon Source



NSF's ChemMatCARS PI/co-PIs and Staff



<https://chemmatcars.uchicago.edu/people/>




Matthew Tirrell, PI and Director, U of Chicago



Mrinal Bera, ASAXS



Wei Bu, Liquid Interfaces



Natalie Chen, User Coordinator



Jason Benedict,
co-PI, U at Buffalo



Ted Betley, co-PI,
Harvard U



Yu-Sheng Chen, Adv.
Crystallography,
Operations Manager



Binhua Lin, Liquid
Interfaces, Exec. Director



Guy Macha, Vacuum
specialist and Tech support



Ka Yee C. Lee, co-PI,
U of Chicago



Mark Schlossman, co-PI,
U of Illinois at Chicago



Mati Meron, X-ray Optics



Charlie Smith, IT




SuYin Wang, Adv. Crystallography

Agenda




| | |
|--|---|
| Overview and Access to NSF's ChemMatCARS | Binhua Lin, <i>NSF's ChemMatCARS</i> |
| Advanced Crystallography Science and Facility | Yu-Sheng Chen, <i>NSF's ChemMatCARS</i> |
| Liquid Interface Scattering Science and Facility | Wei Bu, <i>NSF's ChemMatCARS</i> |
| New Initiative, Anomalous Small Angle X-ray Scattering (ASAXS) Science and Facility | Mrinal Bera, <i>NSF's ChemMatCARS</i> |
| Science Highlights in Advanced Crystallography | Ted Betley, <i>co-PI of NSF's ChemMatCARS, Harvard University</i> |
| Science Highlights in Liquid Interface Scattering | Mark Schlossman, <i>co-PI of NSF's ChemMatCARS, University of Illinois at Chicago</i> |
| Science Highlights in ASAXS | Matthew Tirrell, <i>University of Chicago</i> |

NSF's ChemMatCARS: A Synchrotron X-ray National Facility for Chemistry and Materials Research




Overview and Access




Binhua Lin

NSF's ChemMatCARS
lin@cars.uchicago.edu

NSF's ChemMatCARS – Current Beamline Layout (first experiments in 2002)



One Instrument Operates at a Time



Small-Molecule Adv. Crystallography:

- Charge Density
- Transient States – Photo-crystallography
- Resonant Diffraction
- Structural Dynamics
- Extreme Conditions (High p , Low T)
- Single Micro-Crystal Diffraction
- 3D Δ -Pair Distribution Function

New Initiative

Liquid Interface Scattering:

- Reflectivity
- Grazing-Incidence Diffraction (GID)
- Diffuse Scattering
- Fluorescence near Total Reflection
- Temporal resolution:
1-D Pinhole GID, GIXOS
- GISAXS
- Surface XAS *New Initiative*

Anomalous Small Angle X-ray Scattering (ASAXS):

New Initiative

Beamline Parameters


- Energy range: 5.5 keV to 70 keV
 - X-ray Flux (phot/s/0.1%BW):
 $\sim 4 \times 10^{14}$ @ 10 keV;
 $\sim 1 \times 10^{14}$ @ 30 keV
 - X-ray Energy Resolution: $\Delta E/E \sim 10^{-4}$
- <https://chemmatcars.uchicago.edu/experimental-facility/source-parameters/>

Unusually broad range of techniques serves a wide range of scientific interests

NSF's ChemMatCARS User Access to Beamtime



X-ray beam on 208 days (5000 hours) per year



Open to the scientific community via
APS "General User Proposal" portal



All beamtime proposals peer-reviewed and
ranked by APS external review panels



**Scores determine the allocation of
available beamtime**


NSF's ChemMatCARS Oversubscription Factor is ~ 2
(Oversubscription Factor = Beam days requested over awarded)

NSF's ChemMatCARS Productivity and Broader Impacts



NSF's ChemMatCARS is one of the most productive beamlines at the Advanced Photon Source

Cumulative Growth of User Research Groups Total of ~1400 Users Nationally/Internationally 2004-2019



Publications

<https://chemmatcars.uchicago.edu/publications/>

- 65 -70 Publications annually in peer reviewed journals
- Total Citations 33,772 (as of 9/10/2020)

h-index 91

i10-index 550

Broader Impacts

<https://chemmatcars.uchicago.edu/education-and-outreach/>

- Training the next generation of scientists:
annually 60+% of 200+ users are students/postdocs
- 2 Diversity outreach programs for URM faculty and students
 - Faculty and Student Team Research Award at ChemMatCARS (**FaSTRAC**)
 - Graduate Student Research Award at ChemMatCARS (**GSRAC**)
- Workshops/schools with hands-on experiments and data analysis lectures
3 per year with 40+ participants
- General Public Tours to NSF's ChemMatCARS
~12 per year with ~150 people

New Initiatives at NSF's ChemMatCARS



2-Beamline ChemMatCARS

New techniques at both beamlines – in Adv. Crystallography, Liquid Interfaces, Anomalous SAXS
More capacity to expand in Chemistry and Mat. Research, and extend to Chemical Eng. and Mol. Biosciences


Modernize the **existing beamline** for APS-U (NSF CHE-1834750)

Construct **second beamline** optimized for APS-U (NSF CHE-1836674)

(Two workshops on new initiatives in adv. crystallography and liquid interfaces were held in Nov. 2019)

A Second, Independent Beamline on NSF's ChemMatCARS Sector

2 independent beams to serve 2 concurrent experiments




Apply for Beamtime at NSF's ChemMatCARS




Resources:


- NSF's ChemMatCARS website
- User Coordinator, Natalie Chen,
(nchen@cars.uchicago.edu)
- Lead Beamline Scientists:
 - Adv. Crystallography: Yu-Sheng Chen
(yschen@cars.uchicago.edu)
 - Liquid Interfaces: Wei Bu
(bu@cars.uchicago.edu)
 - ASAXS: Mrinal Bera
(mrinalkb@cars.uchicago.edu)



NSF's ChemMatCARS: A Synchrotron X-ray National Facility for Chemistry and Material Research



Advanced Crystallography Science and Facility




Yu-Sheng Chen
Advanced Crystallography
yschen@cars.uchicago.edu

Advanced Crystallography Science and Facility




Why use **Synchrotron Radiation** for Chemical and Material Crystallography?

- **Intense radiation**
 - **Small crystals $< 2\ \mu\text{m}$**
 - ✓ **Reduces absorption**
 - ✓ **Increases signal to noise ratio**
- **Short wavelength**
 - **High resolution in reciprocal space**
 - **Reduces absorption, anomalous scattering**
- **Monochromatic**
 - **No $k_{\alpha 1}$, $k_{\alpha 2}$ splitting**
 - **Easy data integration**
 - **Better data quality**
- **Tunable energy(Resonant Diffraction)**
 - **Circumvents absorption problems**
 - **Element specific structural information**
- **Fast data collection**
 - **Study of highly air-sensitive samples**



exposure time/120 s
 $d=1.78\ \text{\AA}$




exposure time/0.2 s
 $d=0.69\ \text{\AA}$

Synchrotron data:


- Much less exposure
- Strong data up to **1.78** Å region
- Significant intensity out to **0.69** Å

Advanced Crystallography Science and Facility




Crystal Structure under Extreme Conditions: Beyond Routine Service Crystallography!

•Variable Temperature: 10K ~ 700K




•Variable Pressure: < 10 GPa




Thermal expansion

Compressibility




Phase diagram



Diamond Anvil Cell(DAC)

Advanced Crystallography Science and Facility



Structural Dynamics:

In-situ crystal diffraction

Environmental Control Cell (ECC)


- Vacuum, gas, solution and humidity

Study crystalline materials
under 'real world'
conditions



Photo-Crystallography

- Light-Induced Metastable
State Crystal Structure



Time-Resolved Crystallography




- ☐ Provide an adjustable time window
 - 190 ns and 20 μ s in duration.
- ☐ Easy select single pulses, and longer temporal resolution by adjusting table height.
- ☐ Fully compatible with FPGA

Advanced Crystallography Science and Facility


High energy Pilatus CdTe 1M detector; much better measurements!

3D Diffuse Scattering:

Raw Image(PbTe)




Reciprocal Space




3D- Δ PDF map shows

- ❖ Peaks appear at interatomic vectors
- ❖ Negative-positive-negative profiles, shows correlated displacements




Charge Density:

Deformation Electron Density



Oxalic Acid




- Calculate intermolecular interaction energy.
- Calculate electrostatic potential of the molecule.
- Provide bonding information (topological analysis and source function analysis).
- Characterization of the H-BOND.
- Determine the atomic charge using Atom in Molecular Theory(AIM).
- Extract the "experimental wave function".

Advanced Crystallography Science and Facility

New initiative experiments

- Resonant Diffraction Experiment below 5KeV
- Serial Crystallography for Small Molecule Crystallography



Mail-in Service(remote data collection)



NSF's ChemMatCARS: A Synchrotron X-ray National Facility for Chemistry and Materials Research




Liquid Interface Scattering Science and Facility



Wei Bu
NSF's ChemMatCARS
bu@cars.uchicago.edu

Experimental Hutch




Probe the assembly of molecules, ions, and nanoparticles at liquid interfaces to understand their interfacial distribution, ordering, transport, reactions, and the underlying chemical and physical interactions.

Environmental Chemistry at Liquid Interfaces

- ❖ Separations science (ion extraction) at liquid/liquid interfaces
- ❖ Ion distributions near the interface
- ❖ Surfactant organization for oil recovery and cleanup


(PNAS, 116, 18227, 2019)



Chemistry of Life Processes at Liquid Interfaces

- ❖ Biomolecular recognition
- ❖ Lipid-protein binding
- ❖ Peptide ordering and interactions
- ❖ Interfacial properties of polymers


(PNAS, 111, E1463, 2014)
(Biophysical J., 113, 1505, 2017)



Directed Assembly for Tailored Functionality

- ❖ Functional materials
- ❖ Interfacial complexation for emulsions, films, and coatings
- ❖ Synthesis/ordering of nanoparticles
- ❖ Optical sensor for biomolecule detection
- ❖ Voltage control of interfacial assembly

(PNAS, 115, 9373, 2018)



Why Use X-ray Scattering?

Synchrotron X-ray Scattering is the leading technique for the study of structure at liquid surfaces with high spatial resolution

Molecular organization at the interface


- Molecular ordering and configuration on the sub-nanometer scale
- Structural response to perturbations
pH, electric field, temperature, chemical reactions

Order and disorder within the interface




- In-plane structure of ions, lipids, macromolecules and nanoparticles

Element-specific probes

- Identify elements in complex interfaces
- Monitor variations in interfacial composition
- Local chemical environment at interfaces




Normalized Electron Density



New Beamline Initiatives in Liquid Interfaces

Hybrid interface scattering instrument




<https://chemmatcars.uchicago.edu/experimental-facility/experimental-techniques/liquid-surface-x-ray-scattering/>

New initiatives


Liquid-solid interfaces

(Langmuir 30, 4784, 2014)



- Bio-membrane processes
- Biological cell-surface interactions


Scanning interfacial X-ray probe for heterogeneities at soft interfaces



NSF's ChemMatCARS: A Synchrotron X-ray National Facility for Chemistry and Materials Research



New Initiative,
Anomalous Small Angle X-ray Scattering (ASAXS)
Science and Facility



Mrinal Bera

mrinalkb@cars.uchicago.edu

Scientific Scope


ASAXS: Anomalous Small Angle X-ray Scattering

SAXS done at the energy edges of element of interest

Provides quantitative *element-specific* non-destructive structural characterization in nanometer length scales (1-1000 nm)

How are elements and ions distributed?


Charged Colloids & Macro-ions



[Creative Commons Attribution-Share Alike 3.0 Unported](#)


Ballauf et al, *Annu. Rev. Biophys.*, **40**, 225 (2011)

Nanomaterials




Yarema et al, *Sci. Rep.*, **7**, 1178 (2017)

Polyelectrolytes



Insua et al, *Euro. Polymer J.*, **81**, 198 (2016)

DNA, RNA, & Proteins



Pabit et al, *JACS*, **132**, 16334 (2010)

Expert only technique: No easy to use data reduction and analysis tool


ASAXS Facility @ NSF's ChemMatCARS started in 2017

To develop and bring ASAXS to all users, experts and non-experts alike, in providing them with simple but quantitative element-specific structural characterization tool

ASAXS Facility at NSF's ChemMatCARS

Important features

- Energy Tunable X-rays (5.5-70 keV)
- High X-ray flux $\sim 10^{14}$ at 10keV (ph/sec/0.1%BW)
- Element sensitivity for $Z > 25$
- Software for automated data collection, reduction and ASAXS data modeling



Element sensitivity

| | | | | | | | | | | | | | | | | | | | | | |
|---------------------------|---------------------------|---------------------------|------------------------------|---------------------------|---------------------------|----------------------------|---------------------------|----------------------------|---------------------------|----------------------------|--------------------------|--------------------------|--------------------------|---------------------------|--------------------------|--------------------------|-------------------------|---------------------------|-----------------------|--------------------------|-----------------------|
| H Hydrogen 1.008 | | | | | | | | | | | | | | | | | He Helium 4.003 | | | | |
| Li Lithium 6.941 | Be Beryllium 9.012 | | | | | | | | | | | | | | | B Boron 10.811 | C Carbon 12.011 | N Nitrogen 14.007 | O Oxygen 15.999 | F Fluorine 18.998 | Ne Neon 20.180 |
| Na Sodium 22.990 | Mg Magnesium 24.305 | | | | | | | | | | | | | | | Al Aluminum 26.982 | Si Silicon 28.086 | P Phosphorus 30.974 | S Sulfur 32.065 | Cl Chlorine 35.453 | Ar Argon 39.948 |
| K Potassium 39.098 | Ca Calcium 40.078 | Sc Scandium 44.956 | Ti Titanium 47.883 | V Vanadium 50.942 | Cr Chromium 52.004 | Mn Manganese 54.938 | Fe Iron 55.845 | Co Cobalt 58.933 | Ni Nickel 58.693 | Cu Copper 63.546 | Zn Zinc 65.38 | Ga Gallium 69.723 | Ge Germanium 72.64 | As Arsenic 74.922 | Se Selenium 78.96 | Br Bromine 79.904 | Kr Krypton 83.8 | | | | |
| Rb Rubidium 85.468 | Sr Strontium 87.62 | Y Yttrium 88.906 | Zr Zirconium 91.224 | Nb Niobium 92.906 | Mo Molybdenum 95.94 | Tc Technetium 98.906 | Ru Ruthenium 101.07 | Rh Rhodium 102.91 | Pd Palladium 106.37 | Ag Silver 107.87 | Cd Cadmium 112.41 | In Indium 114.82 | Sn Tin 118.71 | Sb Antimony 121.76 | Te Tellurium 127.6 | I Iodine 126.91 | Xe Xenon 131.29 | | | | |
| Cs Cesium 132.91 | Ba Barium 137.33 | La Lanthanum 138.91 | Hf Hafnium 178.49 | Ta Tantalum 180.95 | W Tungsten 183.84 | Re Rhenium 186.21 | Os Osmium 190.23 | Ir Iridium 192.22 | Pt Platinum 195.08 | Au Gold 196.97 | Hg Mercury 200.59 | Tl Thallium 204.38 | Pb Lead 207.2 | Bi Bismuth 208.98 | Po Polonium 209 | At Astatine 210 | Rn Radon 222 | | | | |
| Fr Francium 223 | Ra Radium 226 | Ac Actinium 227 | Rf Rutherfordium 261 | Db Dubnium 262 | Sg Seaborgium 266 | Bh Bohrium 264 | Hs Hassium 277 | Mt Meitnerium 268 | Ds Darmstadtium 271 | Rg Roentgenium 272 | Cn Copernicium 285 | Nh Nihonium 284 | Fl Flerovium 289 | Mc Moscovium 288 | Lv Livermorium 293 | Ts Tennessine 294 | Og Oganesson 294 | | | | |
| | | | | | | | | | | | | | | | | | | | | | |
| La Lanthanum 138.91 | | Ce Cerium 140.12 | Pr Praseodymium 140.91 | Nd Neodymium 144.24 | Pm Promethium 145 | Sm Samarium 150.36 | Eu Europium 151.96 | Gd Gadolinium 157.25 | Tb Terbium 158.93 | Dy Dysprosium 162.50 | Ho Holmium 164.93 | Er Erbium 167.26 | Tm Thulium 168.93 | Yb Ytterbium 173.05 | Lu Lutetium 174.97 | | | | | | |
| Ac Actinium 227 | | Th Thorium 232.04 | Pa Protactinium 231.04 | U Uranium 238.03 | Np Neptunium 237 | Pu Plutonium 244 | Am Americium 243 | Cm Curium 247 | Bk Berkelium 247 | Cf Californium 251 | Es Einsteinium 252 | Fm Fermium 257 | Md Mendelevium 258 | No Nobelium 259 | Lr Lawrencium 262 | | | | | | |

ASAXS Experiments: Energy dependent SAXS patterns

- Sample shape, size (1-200nm) and size distribution
- Quantitative spatial distribution of elements ($Z > 25$) of choice




ASAXS Data Reduction and Modeling


Tirrell's group (UChicago)

Multilayered
Yb-doped
nanoparticles

side-view



Top-view



ASAXS Measurement

@NSF's ChemMatCARS

